CPR MEDICAL DEVICES INC. Oxylator[®] EMX

The Oxylator® EMX is a patient responsive emergency resuscitation and inhalation management system designed for hospital and EMS personnel. It is simple, safe and effective and delivers optimal levels of oxygen in a resuscitation attempt and dramatically reduces complications associated with time-cycled ventilators and bag valve mask devices. As with other Oxylator® devices, the EMX automatically maintains adequate ventilation for each patient through a patented pressure and flow sensing system that allows for passive exhalation. Unique technology guides the care giver to provide proper airway management.

The Oxylator® EMX

- Reduces the risk of gastric insufflation during mask ventilation due to its low constant flow rate of 30 litres per minute
- Eliminates barotrauma because of a safe pressure limit range of 20 to 45cmH₂O
- Alerts to mask or tube leaks
- Gives the care giver an indication of the patient's lung condition and lets him/her change the ventilatory pattern by manipulating the pressure limits
- Delivers ALL oxygen to the patient; no oxygen is wasted to power the device as in time-cycled devices
- · Frees hands to focus on mask seal to improve airway management
- Informs the care giver of poor lung condition
- Provides manual and continuous cycling modes
- Has an inhalation mode providing "enriched oxygen" for the spontaneously breathing patient
- Assists the shallowly breathing patient
- Eliminates "breath stacking" associated with time-cycled devices
- Offers audible and visual indication of an airway obstruction
- Assures FIO₂ of 1.0 during resuscitation
- Synchronizes automatically with chest compressions
- Reduces the care giver's fatigue
- Always maintains positive airway pressure in continuous cycling mode even at the end of the exhalation phase (2 to 4cmH₂O PEEP)
 - Includes a single use, quick change viral/bacterial filter
 - Is light weight, handy and impact resistant
 - Allows for easy cleaning and assembly
 - Requires minimal training
 - Needs minimal maintenance due to few moving parts
 - •Offers a 5 year warranty
 - · Represents outstanding economic value

Product features and specifications on reverse.

CPR MEDICAL DEVICES INC.

OXYLATOR[®] EMX

The OXYLATOR [®] EMX is sold with five disposable filters, QC report and warranty card. Separate oxygen delivery kits may also be purchased. Please consult your distributor for details.

Feature

Model Weiaht Dimension dia * length Material of housing Required source pressure Required flow (source) **Dead Space** Inspiratory flow rate Minute volume delivered Min. time of oxygen supply I: E ratio PEEP Ventilatory frequency Range of pressure relief Expiratory resistance Suitable body mass range Inhalator flow Oxygen concentration Filter Mask/airway connections Usage temperature range Storage temperature range Obstructed airway warning Oxygen inlet connection Warranty

Specification

OXYLATOR [®] EMX 0.25 kg; 0.55 lbs 2.25 in. * 4.25 in., 57mm * 108mm Acetal between 45 psig to 80 psig, or between 3.0 bar and 5.5 bar min. 30 litres per minute 20 millilitres 30 litres per minute (max) 10 to 12 litres per minute in auto mode Cylinder vol. divided by 12 l/min 1:1 to 1:2 or manually controlled 2 to 4 cm H₂O in auto mode Auto-adjusting to lung capacity in auto mode 20 to 45 cm H₂O , 14.5 to 35 mm Hg approximately 5 cm H₂O in manual mode 10 kg+, 22 lbs + 0 to 15 litres per minute of 100% O₂ 100% during resuscitation mode Single use, disposable 15mm internal / 22mm external - 30 C - + 60 C, - 22 F - +140 F - 40 C - + 70 C, - 40 F - +158 F Rapid cycling, audible and visual DISS, ISO standard or quick connect Five year, maintenance free warranty

SpecificationCons are subject to change.

The OXYLATOR [®] EMX contains patented technology developed and manufactured by:

CPR Medical Devices Inc. 161 Don Park Road Markham, Ontario L3R 1C2 Canada

Tel. 1-416-691-2669 Fax: 1-416-691-7951 E-MAIL: <u>haro@cprmedic.com</u> Website: <u>www.cprmedic.com</u>

U.S. Patent No. 5,230,330 Global patents and further U.S. patents pending